

Ambu *Wireless* Next Generation

Ambu[®] Wireless

Next Generation – discover new independence

Ambu® – Your strong partner

– More than 50 years' experience in medical training

With the innovative development of the first self-inflating resuscitation bag, the Ruben bag, Holger Hesse and Henning Ruben recognized the need for medical demonstration devices. Following the first realistic demonstration training device for using a resuscitation bag, developed by Ambu in 1956, other innovations were developed which enabled the user to learn how to carry out medical procedures safely and with long-term success. With the launch of the Phantom, Ambu entered the field of professional medical training in 1959 and has kept its finger on the pulse ever since, continuously developing and improving its product range. It is important to Ambu that the products are geared perfectly to the needs of the user – that's why they are developed in close cooperation with doctors, nurses, paramedics and medical research groups.



All Ambu manikins are continuously adapted and enhanced to meet the changing requirements of modern professional training. This means taking into account not only the latest medical techniques, but also the change in training concepts and teaching methods. In addition to the anatomical and physiological authenticity of the manikins, quality and durability above all lie at the heart of these developments. The functionality and durability of all Ambu manikins is tested extensively. The mechanical construction of the manikins also keeps maintenance costs down. Together with good customer service and comprehensive spare parts availability, Ambu manikins offer a long-term, cost-effective and reliable solution for professional medical training, as well as basic first-aid training.



1959
Ambu
Phantom



1977
CPR Pal



1985
Simulator II



1991
i.v. Trainer



1993
Defib Trainer



1973
Simulator



1982
Intubation
Trainer



1987
AmbuMan
CPR Software



1992
Baby



1994
MultiMan

The future begins with W ...

... for wireless. Ambu showed this back in 2009 by being the first provider in the manikin market to develop a fully integrated wireless solution. Ambu Wireless Technology was born and now allows a new level of scenario-based training.

In 2016, Ambu went one step further. Based on customer feedback, we have enhanced Ambu Wireless Technology and created a completely new concept that is perfectly tailored to the needs of modern course concepts for Advanced Life Support (ALS).



* 1996
MegaCode
Trainer

* 1999
Mini MegaCode
Trainer

* 2003
M MegaCode

* 2009
Ambu
Wireless

* 2015
AmbuMan Compression
AmbuMan Basic
Kit i.o. Access

* 1997
Airway Management
Trainer

* 2001
UniMan

* 2008
Airway Man

* 2011
AmbuSam

* 2016

Ambu Wireless
Next
Generation

2016

We're turning everything on its head:

From WWW to MMM

Ambu® Manikin Management Module

At the heart of the concept is the AmbuMan, which features the latest integrated Ambu Wireless Technology.

Control, documentation and monitoring of the manikin takes place wirelessly via the browser-based Ambu Manikin Management Module.

The Ambu Manikin Management Module is the central interface where all of the data and information converge: compression rate, compression depth, hand position, ventilation volume and gastric insufflation.

Infinite

According to the 2015 AHA Guidelines, realistic scenario-based training is essential – unleash your creativity and discover the new possibilities the AmbuMan Advanced has to offer!




Realistic

Familiar, realistic anatomy. A highlight of the AmbuMan Advanced is the representation of a patient monitor with 12 lead ECG display.

Easy

Open the Ambu Manikin Management Module in the browser – without any additional software installation!



Flexible

Connect from virtually any PC, tablet or smartphone – wirelessly and platform-independent!



Secure

Complete separation of the network configuration level and the training/control level – ready to start in minutes!

Unlimited

Parallel multiple access of different terminal devices for performance visualization is possible. Furthermore, the monitoring module can be adapted to your individual requirements.

Individual

The option to provide individual evaluation in group training ensures effective training. Each participant receives their own personal results (up to four participants).

Independent

The AmbuMan Wireless features a license-free self-training module as standard which allows training stations to be set up and non-trainer-specific training to be conducted at any time.



Ambu® Wireless Manikins

– Realistic training, more flexible than ever before

Ambu Man Wireless

The all-rounder for CPR: can be used with both healthcare professionals and for training lay rescuers. Features the license-free self-training module and Ambu hygienic system.



Ambu Man Defib

The specialist in early defibrillation. Offers the basic functions of the AmbuMan Wireless with additional defibrillation option.



Ambu Man Airway Wireless

The specialist in airway management: for training in both endotracheal intubation and alternative techniques like laryngeal mask or laryngeal tube.

Ambu Man Advanced

The professional for professionals. Ambu's new generation of manikins, perfectly adapted to the training of Advanced Life Support (ALS).




Ambu® Wireless Next Generation

– Discover NEW independence

The network configuration level has been separated completely from the training and control level. This provides added convenience and security for your training. Once the manikin has been set up, it is ready to start within minutes. Complete separation makes unintentional access at network

level during training impossible, and unwanted breaks in training are no longer possible.

At the training and control level, the different modules can be configured individually and adapted to the requirements of your course concepts.



Flexible and intuitive control

Discover the new flexibility and platform independence. Ambu Wireless supports all modern network connection options that the TCP/IP network offers. You can connect with your Ambu manikin wirelessly from virtually any PC, tablet or smartphone. Furthermore, all AmbuMan Wireless Manikins are fitted with NFC technology, which allows a fast and uncomplicated connection.*

The wireless star – connect your manikins


To conduct a group training session, you connect several manikins together to form what is known as a star network: all of the manikins are connected to a central hub, the host manikin. For training, we recommend connecting five manikins to one host to ensure that the displays remain clear.



WLAN



LAN



*Requirement: your terminal device must support NFC.

Ambu® Wireless

– Discover infinite scenarios

With Ambu Wireless you can provide targeted training. Create your own scenario down to the finest detail. Simulate resuscitation in conjunction with standard emergency procedures, like transportation and repositioning. Even airway management in a trapped car can be realistically simulated and practiced. Meanwhile, you can supervise the scenario from a safe distance. Thanks to the wireless transmission of the most important parameters of the CPR, the training can be conducted with a high level of mobility.



Wireless mobility and flexibility

As the trainer, you can observe the training session from the sidelines, without having to pay attention to wiring. The manikin is the centerpiece of the training, and after an initial one-time configuration in the network configuration module, the manikin is ready for the next training session within minutes. You can concentrate fully on the training and unwanted access to the network becomes a thing of the past.

If training is taking place in a building, you can use the training rooms according to your requirements and adapt the network to the space available; or transmit the training parameters to other rooms onto different terminal devices simultaneously, for example onto a Smart Wall or different tablets. If you are conducting training outdoors, Ambu Wireless offers the highest level of mobility, with a range of up to 50 meters.

CPR feedback in real time and debriefing

The Ambu Manikin Management Module allows a new level of debriefing. The option to provide individual evaluation in group training ensures effective training in all course content. Each participant receives their own individual results.

With the Ambu Manikin Management Module, you receive immediate feedback on the training

parameters. For example, you can display the compression depth and ventilation volume on the stomach instrument and in the monitoring module of the Ambu Manikin Management Module. You can also incorporate the realistic patient monitor into your training, which gives the trainee feedback on the physiological parameters and a 12 lead ECG display.

Durability and robustness

The Ambu training product line offers innovation based on tried-and-tested quality. Ambu manikins are highly regarded worldwide for their robust design and the high quality of materials used in their construction.

Your creativity

With all of the ideas you have for creating realistic training, scenario-based training comes alive. The days of forming a circle with a manikin in the middle are long gone. Discover new possibilities.

The 2015 AHA Guidelines outline the curricular content and requirements of resuscitation training. As a manufacturer of professional manikins, we show you the possibilities of implementing this content in your training sequences using our products.

Simulation and training of realistic scenarios

Simulating situations as realistically as possible is seen as an essential part of resuscitation training according to the AHA core educational components. We share this view; with the wireless generation, you can create realistic scenarios and better prepare your participants for a genuine emergency. There are virtually no limits to your creativity or the individual configuration of emergency settings.

Debriefing – to improve CPR

Debriefing following training, during which objective facts are discussed openly, can help improve CPR. With the Ambu Manikin Management Module, you can start an in-depth discussion by referring to all of the events that have been recorded. Illustrate to the participants which measures were successful and which activities or situations that need improvement.

Increased focus on non-technical skills (NTS)

Using scenario-based training, important skills like leadership, teamwork, task allocation and structured communication can be taught creatively and comprehensively.

Refresher courses to maintain skills

After a while, knowledge and skills acquired through training begin to decline. The Ambu Manikin Management Module offers a database function. This allows you to save your participants' individual performance and compare them with each other in subsequent training sessions.

Effective preparation for training

Alongside your individual preparation for the training, as the instructor you can also structure and allocate groups of participants before the actual training day using the course planner integrated into the Ambu Manikin Management Module.



Ambu Manikin Module – for evaluating and debriefing

"The provision of feedback and/or debriefing is a critical component of experiential learning. Feedback and debriefing after skills practice and simulations allow learners the opportunity to reflect on their performance and to receive structured feedback on how to improve their performance in the future."¹

AmbuMan Advanced for ALS training

"Standard manikins continue to be an appropriate choice for organizations that do not have this capacity."¹

Ambu® Manikin Management Module

– The "software" that can do MORE



A new level of debriefing

For the first time, the Ambu Manikin Management Module makes it possible to carry out individual appraisals in group training. Using the switch function, every trainee is assigned their individual training content. This means that every trainee receives their individual appraisal and results (up to four participants).



Customizable ECG simulation

The Ambu Manikin Module contains an ECG module which you can use to simulate an ECG, cardiac arrhythmias and artifacts wirelessly on the AmbuMan Advanced and AmbuMan Defib. It is particularly important that the ECGs are physiologically accurate and that they will be correctly detected by an external defibrillator. The module offers extended options and you can change existing ECG rhythms or create your own.



Interaction

All parameters can be exported/imported. This gives you the option of changing the parameters between different Ambu Wireless manikins and provides a way to interact between different centers.

AMBU W
Man
Manag
Mod
NEXT
GEN

WIRELESS GENERATION

Ambu
Manikin
Management
Module



Event management

Select the events from the algorithms, such as "Speak to and shake" or "Tilt head and lift chin", to fix them within the time-frame of the resuscitation. During the resuscitation and as part of the evaluation or subsequent review, you can measure and discuss the times between the events. Use the software to create your own events, like "Defi 150J" or "Ami 300 mg".




Playback of sound files

On the AmbuMan Advanced it is possible to play back sound files, such as breathing in the scenario. You can also install your own individual sound files onto your manikin.



Customizable displays

Make the most of the flexible display options within the monitoring module. The display can be adjusted flexibly to the requirements of the audience.





Patient monitor

– Realistic representation of physiological parameters

A highlight of the Ambu® Man Advanced is the representation of a realistic patient monitor with the following display options: **12 lead ECG, heart rate, blood pressure, oxygen saturation, capnography and capnometry.** The physiological

parameters can be configured individually by the trainer and changed during the training session, and are immediately visible to the trainee.



On your mark – get set – go!




– Self-training module

Nothing is more important than revising learned content and practicing the basics of CPR until they become automatic and routine in emergency situations.

The Ambu® Man Wireless, Ambu® Man Defib and Ambu® Man Airway Wireless feature a self-training module as standard which allows training stations to be set up and non-trainer-specific training

to be conducted at any time. This means that resuscitation teams can train continuously and trainees can prepare for their tests.

This module comes as standard and is license-free with the AmbuMan Wireless, AmbuMan Defib and AmbuMan Airway Wireless, allowing cost-effective and continuous self-training.



Leistung	Bewertung	Gesamt	Durchschnittswert	Empfehlung	Evaluierung
Beatmungsvolumen pro Minute	●		2.6 L/Min	1.5 - 3.6 L/Min	
Herzzeitvolumen	●		1.8 L/Min	2 - 4.8 L/Min	
Beatmungsrate	●		3 /Min	3 - 9 /Min	
Beatmungsvolumen	●		650 mL	500 - 600 mL	
Magenblähung	●	0	0		
Kompressionsrate	●		62 /Min	100 - 120 /Min	
Kompressionstiefe	●		55.6 mm	50 - 60 mm	
Verhältnis Kompression/Entspannung	●		63/37	50/50	
Falsche Handposition	●	0	0		
Fehlende Entlastung	●	0	0		
Verhältnis Kompression/Atmung	●		62/2	30/2	

Ambu[®] Man Wireless Manikins

Specifications

Mouth-to-mouth ventilation

Mouth-to-nose ventilation

Bag valve mask ventilation

Use of Guedel airways

Airway management: endotracheal tube, laryngeal tube, laryngeal mask, other SGAs

Jaw-thrust maneuver

Gastric insufflation

Ambu hygienic system

Lung bag

Face masks

Chest compression

Hand position display

Adjustable chest stiffness

Pulse

Placement of defibrillation electrodes

Defibrillation (up to 360 joules)

ECG display

Transthoracic pacing option

i.v. access

i.o. access

Advanced feedback display (ventilation volume, compression depth, hand position, gastric insufflation)

Ambu Wireless Technology

Ambu Manikin Management Module

Platform-independent access and control options

Browser-based access and control option

Parallel control of up to 6 manikins

Scenario-based training (individually programmable)

Manual event management

User-defined programming of ECG rhythms

Debriefing (individual documentation of up to 4 participants)

Playback of sound files e.g. breathing

Parameter export/import

Patient monitor display of physiological parameters like 12 lead ECG, heart rate, blood pressure, oxygen saturation, capnography and capnometry.

Performance display (customizable)

Parallel display of up to 6 manikins

Parallel multiple access of different terminal devices

License-free self-training

Secure separation of network level and training/control level

Ambu Manikin Management Module

Control module

Monitoring module

Self-training module

Network configuration module



Targeted training

Easy

Infinite

Security

Flexible

Independence

wireless

Convenience

Platform independence

Realistic

Interactive

Individual

Efficient