
Instruction for Use


Flow Splitter


English	4-5
Español	6-7

The Flow Splitter is designed to uniformly split the Fuser Pain Pump and/or the Block Pain Pump flow and to enable the simultaneous administration of medication through two catheters.

1. Connect the administration set without the catheter and fill the pump.*
2. (a) Remove the caps from the Flow Splitter

- (b) Connect the Flow Splitter to the pump exit port (luer lock) and connect both catheters to the Splitter's two exit ports (luer locks).

NOTE: Verify that the luer lock connectors are tightened.

3. Perform priming to the administration set together with the Flow Splitter and the two catheters according to the IFU.*
4. Set the pump flow rate; consider the below flow-rate table and the IFU.*

Regulator Flow Rate	Catheter's Flow Rate
3 ml/hr	1.5 +1.5 ml/hr
4 ml/hr	2.0 + 2.0 ml/hr
6 ml/hr	3.0 +3.0 ml/hr
8 ml/hr	4.0 +4.0 ml/hr
10 ml/hr	5.0 +5.0 ml/hr
12 ml/hr	6.0 +6.0 ml/hr
14 ml /hr	7.0 +7.0 ml/hr

* Refer to Physician Guidelines and Instructions included with the pump set.

NOTE

- **Single Use product, do not reuse**
- **The Flow Splitter is for single use only. Do not re-sterilize.**
 - (a) re-use might result with patient contamination.
 - (b) re-sterilization will result in poor device performances that do not meet the declared tolerances.
- **Sterile, unless sterile packaging is broken. Do not use the product if any damage have occurred to the sterile packaging.**
- **The product is sterile by Ethylene Oxide - STERILE EO.**

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