

# AmbuMan<sup>®</sup> School

## DATASHEET

AmbuMan School is a training manikin dedicated to teaching basic life support techniques. The manikin can be used to realistically train on external chest compressions and place AED pads.

The unique hygienic system prevents cross-contamination when demonstrating and training the principals of ventilation through the mouth and nose.

The AmbuMan School is intended to be used to educate students in a school setting under the supervision of an instructor/teacher.

### Conforms to guidelines

The AmbuMan School allows instructors to perform CPR training following the latest CPR guidelines.

### Compact and portable

Lightweight and compact, the AmbuMan School is easy to transport and set up.

### Easy AED training

Due to the closed shaped torso, defibrillator pads can be placed and simulated AED training can be performed.



## KEY BENEFITS

Ambu hygienic system to prevent cross-contamination

Compact and portable training manikin

Easy training of chest compressions

Closed torso

Easy placement of AED pads

10 manikins for optimal class training

# Ambu

## SPECIFICATIONS

Weight	0.8 kg
Length	45 cm
Max. compression depth	6 cm

## MATERIALS

Product	Material
Straps	Nylon
Base plate	Polycarbonate
Bellow	Thermoplastic Elastomer
Chest	Polypropylene
Head	Silicone
Face piece	Polyvinylchloride
Lung bag	Polyvinylchloride
Head bag	Polyethylene

## ORDER INFORMATION


Item No.	Description
A268 007 000	AmbuMan School (Set of 10)
268 000 571	Face piece AmbuMan School (Set of 10)
268 000 572	Head bag AmbuMan School (Set of 100)
268 000 573	Accessories AmbuMan School (Set of 10)



Closed torso



Placement of AED pads



Ambu hygienic system



Accessories

**Ambu** FOREVER FORWARD

Ambu, Inc.  
6721 Columbia Gateway Drive  
Suite 200  
Columbia, MD 21046  
Tel. 800 262 8462  
Fax 800 262 8673  
ambuUSA.com

Questions? Contact your Ambu representative at  
**800-262-8462** or visit **www.ambuUSA.com**

CE