

**THINK FAST
THINK CLEAN
THINK SINGLE-USE**



Ambu® aScope™ Gastro

Ambu

KEY BENEFITS

- Available at all times and in any setting
- A new sterile gastroscope for every patient
- Consistent quality, feel and performance
- No hidden costs

TECHNICAL SPECIFICATIONS

Optical System

Field of View	140°
Depth of Field	3 - 100 mm
Illumination Method	Dual LEDs
Image Enhancement	Yes (Advanced Red Contrast)

Insertion Tube

Distal Tip Outer Diameter	9.9 mm
Working Channel Inner Diameter	2.8 mm
Working Length	1030 mm

Bending Section

Angulation (Up/Down/Left/Right)	210°/90°/100°/100°
---------------------------------	--------------------

Endoscope Channels

Working Channel Inner Diameter	2.8 mm
Auxiliary Water Jet Channel	Yes

Weight

Weight	650 g
--------	-------

Sterilization


Sterilization	ETO
---------------	-----

Storage and Transportation

Transportation Temperature	-10 - 55° C
Storage Temperature	10 - 25° C
Relative Humidity	10 - 95 %
Atmospheric Pressure	50 - 106 kPa

ORDERING SPECIFICATIONS


Description	Reference number	Units/Box
aScope Gastro	483001000US	6 per box



210° Retroflexion of the aScope Gastro



The Compatible Ambu aBox 2 Display and Processing Unit



Ambu FOREVER FORWARD

Ambu Inc.
6230 Old Dobbin Lane
Suite 250
Columbia, MD 21045
Tel. 800 262 8462
Fax 800 262 8673
www.ambuUSA.com

Questions? Contact your Ambu representative at **800-262-8462** or visit www.ambuUSA.com/single-use-gastroscope

US Rx Only