



## The Standard in Resuscitation Bags.

### For Manual Ventilation of Adult, Children, and Neonates

The Ambu SPUR II (Single Patient Use Resuscitator) is designed for manual ventilation of neonates through to adults.



#### Key FEATURES

- Unique single-shutter valve system for reliable functionality
- Integrated handle for user comfort and uniform compression
- Thin-walled compression bag allows for lung compliance and “feel”
- SafeGrip™ surface for secure handling in stressful environments
- Medication port permits quick medication delivery without disconnecting the Ambu SPUR II from the ET tube
- Easy attachment of manometer and PEEP valve

## Ambu SPUR II Accessories



### Ambu Peep Valve



#### Key FEATURES

- Adjustable between 1.5-20 H<sub>2</sub>O cm
- 30 mm connector and fits all sizes of the Ambu SPUR II, as well as most competitive resuscitators

### Ambu CO<sub>2</sub> Detector



#### Key FEATURES

- Designed for fast, easy airway assessment
- Verifies proper endotracheal tube placement
- Helps prevent unwarranted intubations
- Confirms cardiac output during CPR
- Expansive 2" diameter
- 15x22 mm allows the Ambu CO<sub>2</sub> Detector to be attached to any brand resuscitator

### Ambu Manometer



#### Key FEATURES

- Unique Barrel design offers 360° viewing capabilities
- Monitors pressure between 0-60 cm H<sub>2</sub>O and fits all resuscitators equipped with a manometer port
- The slim design makes it easy to monitor pressure with minimum obstruction of the of the patient's face
- Single patient use
- All materials are latex and PVC free

### Filter Media



#### Key FEATURES

- Various models offered:
  - Inspiratory & Expiratory
  - HEPA & Bacterial/Viral
- Fitting Dimensions:
  - 30 mm x 30 mm
  - 15 mm x 22 mm

**Ambu Inc.**  
6740 Baymeadow Drive  
Glen Burnie  
Maryland 21060  
Phone:(800) 262-8462  
Fax: (800) 262-8673  
[www.ambu.com](http://www.ambu.com)

RESPIRATORY CARE

