

SIMPLIFY YOUR SWALLOW ASSESSMENTS

FEES exams with the
Ambu® aScope™ 4 RhinoLaryngo



Ambu



SIMPLIFY YOUR FEES EXAMS

The single-use aScope™ 4 RhinoLaryngo Slim is always ready and available when you need to do a Flexible Endoscopic Evaluation of Swallowing (FEES) exam. Unlike reusable laryngoscopes—where you have to consider the volume and location of patients and plan accordingly—single-use scopes enable you to perform exams on patients where they are and when they need it.

aScope 4 RhinoLaryngo Slim and aView 2 Advance display system is a lightweight, portable solution you can easily transport into ICUs, hospital wards and acute care facilities, and is simple to position near the patient.

Easy to set up for use - you simply remove the endoscope from its sterile packaging, connect it to the full-HD aView™ 2 Advance monitor, and you see a live image within seconds. Additionally, the aView system can run on battery power for over two continuous hours - enabling clinicians the flexibility of exams and documentation in virtually any setting.

A NEW VIDEO ENDOSCOPY SOLUTION: Low capital investment, no repairs, and distal chip quality every time

Along with low capital investment compared to reusable scopes, the single-use aScope™ 4 RhinoLaryngo Slim also eliminates hidden costs associated with reusables, such as cleaning supplies, maintenance, repairs, and long-term contracts. By investing in single-use - expenses are predictable and ideal for expanding procedure volumes and points of care without more scopes and repairs.

The aScope 4 Rhinolaryngo was found to be more cost effective per use than reusable flexible laryngoscopes based on a five-year life span.*

*R. Walczak, M. Arnold, J. Grewal, X. Yuan, A. Suryadevara, and H. Marzouk, "Reusable vs disposable nasopharyngolaryngoscopes : Cost analysis and resident survey," Laryngoscope Investig. Otolaryngol., no. September, pp. 1-6, 2020, doi: 10.1002/liv.2.500.



CLEAR AND EASY DOCUMENTATION

The aScope™ 4 RhinoLaryngo Slim solution enables you to easily record and store images and videos from an assessment. These are available for later viewing, including frame-by-frame playback, which you can review with patients and peers, or use to support post-assessment documentation. Additionally, the aView 2 Advance platform has outputs that allow for connection to swallowing/FEES-centric documentation system software.

Finally, the aView 2 Advance can connect with internal PACS/EHR systems for communication and data transfer.



Here's how the aScope 4 RhinoLaryngo Slim can make a difference in FEES procedures:

- Provides immediate access, cuts down on patient wait time, and is ideal on weekends when skeleton staffing can mean limited barium swallow test availability
- Eliminates reprocessing following exams, along with point-of-care difficulties and endoscope transport
- Is lightweight and improves ergonomics
- Cuts out capital or service/repair expenses after repeated usage of traditional endoscopes
- Moves away from fiber-optic/camera setups and bulky, stationary video towers
- Enables frame-by-frame video replay of swallow physiology
- Offers objective data that does away with over-diagnosis and under-diagnosis of aspiration risk
- Offers highly portable, compact HD display that is plug and play, mounts easily to IV poles, and operates on battery power when needed
- Ambu has a nation-wide network of Clinical Training Specialists to support new users in training and hospital/clinic installation, including on-site education and training.

**For more information,
please visit ambuUSA.com/FEES**

Ambu, Inc.
6721 Columbia Gateway Drive, Suite 200
Columbia, Maryland 21046
Tel. 800 262 8462
Fax 800 262 8673
ambuUSA.com

US: Rx only