

UNLEASH THE POWER OF SINGLE-USE

About Ambu

Since 1937, Ambu has been rethinking medical solutions to save lives and improve patient care. Millions of patients and healthcare professionals worldwide depend on the efficiency, safety and performance of our single-use endoscopy, anesthesia, and patient monitoring solutions.



Ambu® Product Portfolio

Ambu

ANESTHESIOLOGY & RESPIRATORY

Our latest innovations for anesthesiologists will help save lives and improve patient care. This full portfolio of quality anesthesia products crosses over a variety of applications, from everyday usage to exceptional situations.



Masks

The King Mask is Ambu's premium single-use disposable anesthetic face mask. The anatomically shaped mask features a thin cushion that is tacky to the touch which supports a tight seal. The top of the cushion is ribbed for a no-slip grip for the clinician. All masks are DEHP (phthalate) free, latex free, and MRI safe. Pediatric masks are scented.



Breathing Circuits

The innovative designs of King Single Limb circuits improve patient care. The Unique F2 filter protects both patient and machine (inspiratory and expiratory ports) with a single, dual-chamber, and filter. Provides bacterial and viral filtration of 99.99% or greater. The ability to shape and expand flex tubing to ideal length reduces torque on the airway device. Coaxial design helps keep patient normal thermic. Flex2 reduces hospital waste.



Breathing Bags

Ambu's King Breathing Bags provide a secure attachment to the machine as well as an excellent feel. The contour hourglass shaped bag allows for easier hand use. Our colored bags help distract pediatric patients and are latex free.



Ambu® SPUR® II Resuscitators

Ever since the invention of the world's first resuscitator in 1956 – the now iconic Ambu bag – Ambu has become the epitome of reliability and convenience.



Laryngeal Tubes

The Ambu® King LTS-D™ is a dual cuff alternative airway device that includes a separate gastric channel to support stomach decompression and suctioning. Reliable and easy to use, the King LTS-D is a powerful tool for improving out-of-hospital cardiac arrest outcomes in adults.*

* Life-saving-related claims are out of the FDA product code scope. This is a scientific communication based on this specific clinical study. It is not cleared by FDA SIUU communication guidance.



Laryngeal Masks

The AuraGain is Ambu's 2nd-generation laryngeal mask. It satisfies three fundamental airway management needs by integrating gastric access and intubation capability in an anatomically curved single-use device that facilitates rapid establishment of a patent airway.

REVOLUTIONIZING RESUSCITATION

Ambu debuted the world's first self-inflating resuscitator, designed for manual ventilation, in 1956. Since then, the Ambu Bag has become synonymous with saving lives.



PATIENT MONITORING

Our patient monitoring and diagnostics product line includes innovative solutions for your cardiology, neurology, and intraoperative monitoring needs.



ASK US ABOUT
OUR SINGLE PATIENT
PACKAGING OPTIONS
FOR CUP ELECTRODES

Neurology


Ambu has worked for decades with neurophysiology experts to co-develop tools that meet the quality you demand.



Scan to download
the Neurology
Catalog




EEG Cup Electrodes



EMG Guided Injections




Surface Electrodes




EMG Needles

Intraoperative Monitoring

Protect patients during surgery by monitoring their nervous system with our line of surface electrodes, probes, and subdermal needles.



Subdermal Needles



Twisted Pair



Probes



Surface Electrodes

Cardiology

Since 1976, when the first single-use Ambu® BlueSensor™ ECG electrodes were introduced, we have delivered quality solutions matched to your specific needs. Our full portfolio, with its wide range of sizes, shapes, gel types and backing materials, provides devices for every situation.



Scan to download
the Cardiology
Catalog



Neonatal



Pediatric



Adult






AIRWAY MANAGEMENT

The inventor of the world's first self-inflating resuscitation bag offers a full portfolio of products for airway management professionals to help save lives and improve patient care. Explore the clinical evidence and learn more about our products and resources.



VivaSight™ 2 DLT

VivaSight 2 DLT is the World's first single-use Double Lumen Tube with integrated camera allowing for continuous visualization during One Lung Ventilation procedures.



Fuji UNIBLOCKER™

UNIBLOCKER is the only bronchial blocker that incorporates a steel mesh polyurethane-covered shaft for direct control of the blocker into the targeted bronchus. This control allows for simple and intuitive placement of the blocker.



AMBU'S COMMITMENT TO REDUCE OUR ENVIRONMENTAL FOOTPRINT

Our Mission — To create solutions that improve health-care for patients and professionals with the smallest environmental footprint possible.

- **Integrated circular approach** from product design to manufacturing and device disposal
- **Renewable materials** that reduce impact of plastic in products and material for packaging
- **Take-back and recycling processes** to keep products from landfills and recycle material for other products

Recyclable Packaging



All our secondary packaging is
***100% recyclable**



Rethinking packaging to **utilize renewable or recycled materials**



Initiated **implementation of bioplastics** in protective parts

* Secondary packaging includes protection tubes but not Tyvek pouches

Carbon Emissions & Climate Change



Near-term carbon reduction targets **validated by Science Based Targets Initiative (SBTi)** in Dec. 2023



Committed to **reducing Scope 1 and 2 greenhouse gas emissions by 75%** by FY 2029/30



Working closely with key suppliers to **reduce Scope 3 emissions** in our value chain



Target goal of **82% of suppliers setting science-based targets by FY 2026/27**

THE AMBU ECOSYSTEM

We provide healthcare supply chain executives with a comprehensive solution to improve patient safety, operational efficiency, and cost-effectiveness. Our rapid one-step setup with an ultra-mobile system offers ease of use across multiple hospital settings, from the outpatient clinic to inpatient cases.

aCart™ Compact

Ambu's new aCart™ Compact is a modern, mobile workstation providing easy storage for our single-use endoscopy portfolio in a small footprint – convenient, portable, versatile and hygienic.



aCart™ Plus


aCart™ Plus offers a stable and convenient workstation for healthcare professionals who perform procedures in the bronchoscopy suite or other endoscopy suites.



aScope™ 5
Cysto HD



aScope™ 5 Uretero



aScope™
Duodeno 2



aScope™
• Gastro
• Gastro Large




By 2025, all Ambu single-use endoscope handles will be made from bioplastic.



aScope™ 5

- Broncho HD 5.6/2.8
- Broncho HD 5.6/2.8 Sampler Set
- Broncho HD 5.0/2.2
- Broncho HD 5.0/2.2 Sampler Set
- Broncho 4.2/2.2
- Broncho 4.2/2.2 Sampler Set
- Broncho 2.7/1.2



**Ambu®
VivaSight 2 DLT**



aScope™ 4

- Broncho Slim
- Broncho Regular
- Broncho Large
- Broncho Sampler Set



aScope™ 4 Cysto



aScope™ 4

- RhinoLaryngo Slim
- RhinoLaryngo Intervention



**aView™ 2
Advance**

aBox™ 2

SINGLE-USE ENDOSCOPY

Ambu created the world's first single-use flexible endoscope in 2009. We now leverage more than a decade of single-use flexible endoscopy experience in multiple clinical areas, with more to come in the near future. Join the world leader in single-use endoscopy on this journey.



Pulmonology

With our single-use bronchoscopy solutions, you always have access to a new sterile bronchoscope with consistent quality, for basic and advanced procedures.

aScope™ 4 Broncho Slim, Regular, Large
aScope™ 4 Broncho Sampler Set
aScope™ 5 Broncho HD 5.6/2.8, 5.0/2.2
aScope™ 5 Broncho HD 5.6/2.8, 5.0/2.2 Sampler Set
aScope™ 5 Broncho 4.2/2.2, 2.7/1.2
aScope™ 5 Broncho 4.2/2.2 Sampler Set
VivaSight™ 2 DLT



Gastroenterology

The Ambu GI portfolio of endoscopes offers an ideal solution for critical care and surgical endoscopy. With unmatched efficiency, the right scope is always available exactly where and when needed, empowering your team to deliver exceptional patient care without delay.

aScope™ Gastro
aScope™ Gastro Gastro Large
aScope™ Duodeno 2



ENT (Ear, Nose, Throat)

Ambu offers single-use rhinolaryngoscopes with high-quality imaging, multiple connectivity options for the outpatient clinic and easy portability for when you're called to a consult.

aScope™ 4 RhinoLaryngo Slim

aScope™ 4 RhinoLaryngo Intervention



Urology

Ambu's single-use urology solutions increase your productivity and simplify your workflow.

aScope™ 4 Cysto

aScope™ 5 Cysto HD

aScope™ 5 Uretero



ORDERING SPECIFICATIONS

Item No.	Description
476101000	aScope™ 4 Broncho Slim 3.8/1.2
476001000 (Canada)	aScope™ 4 Broncho Slim 3.8/1.2
477101000	aScope™ 4 Broncho Regular 5.0/2.2
477001000 (Canada)	aScope™ 4 Broncho Regular 5.0/2.2
478101000	aScope™ 4 Broncho Large 5.8/2.8
478001000 (Canada)	aScope™ 4 Broncho Large 5.8/2.8
477501000	aScope™ 4 Broncho Regular Sampler Set
478501000	aScope™ 4 Broncho Large Sampler Set
624001000US	aScope™ 5 Broncho 2.7/1.2
620001000US	aScope™ 5 Broncho 4.2/2.2
621001000US	aScope™ 5 Broncho HD 5.0/2.2
622001000US	aScope™ 5 Broncho HD 5.6/2.8
621002000US	aScope™ 5 Broncho Sampler Set 5.0/2.2
622002000US	aScope™ 5 Broncho Sampler Set 5.6/2.8
601101000	aScope™ 4 Cysto, Standard Deflection
601001000 (Canada)	aScope™ 4 Cysto, Standard Deflection
600101000	aScope™ 4 Cysto, Reverse Deflection
600001000 (Canada)	aScope™ 4 Cysto, Reverse Deflection
603101000	aScope™ 5 Cysto HD, Standard Deflection
602101000	aScope™ 5 Cysto HD, Reverse Deflection
604101000	aScope™ 5 Uretero, Standard Deflection
605101000	aScope™ 5 Uretero, Reverse Deflection
510101000	aScope™ 4 RhinoLaryngo Slim
510001000 (Canada)	aScope™ 4 RhinoLaryngo Slim
512101000	aScope™ 4 RhinoLaryngo Intervention
512001000 (Canada)	aScope™ 4 RhinoLaryngo Intervention
486104000	aScope™ Duodeno 2
483001000US	aScope™ Gastroscope
582101000	aScope™ Gastro Large
412351000	Ambu® VivaSight 2 DLT 35 Fr Left
412371000	Ambu® VivaSight 2 DLT 37 Fr Left
412391000	Ambu® VivaSight 2 DLT 39 Fr Left
412411000	Ambu® VivaSight 2 DLT 41 Fr Left
505001000US2	aBox™ 2
405011000US2	aView™ 2 Advance (Gen. 2) monitor
KU.9713.800	aCart™ Compact
KU.9740.800	aCart™ Plus

Call your Ambu customer service representative at 800-262-8462 for ordering options on the following products

Cardiology
Neurology
Circuits & Masks
Resuscitators
Airway
EMS

Ambu FOREVER FORWARD

Ambu, Inc.
6721 Columbia Gateway Drive
Suite 200
Columbia, MD 21046
ambuUSA.com

US Rx only

PUB - 001585 - V05 - 11/2024 - Ambu USA. Technical data may be changed without notice.