


Ambu® Head Wedge®

DATASHEET

The unique Ambu Head Wedge is designed to immobilize the patient's head safely and effectively when placed on a backboard.

Its innovative design provides a secure fit on the patient in four quick and easy steps. This saves time when every second counts.



KEY BENEFITS

Waterproof and durable

Radiolucent to safely immobilize patient during x-ray

Stores flat to reduce valuable storage space

May be stored attached to the back board

Fully adjustable to fit all sizes from pediatrics to adults

Fits all backboards


Strong, self-adhesive head strap ensures a secure hold without any discomfort to the patient

Ambu

SPECIFICATIONS

Description


Accommodates head sizes	50 cm to 63 cm in circumference (19.7" to 24.8")
Dimensions	(L x W x T) 400 x 165 x 12 mm (15.8" x 6.5" x 0.5")
Weight	100 g (3.5 ounces)
Recommended operating temperatures	-5°C – +40°C (-23°F – + 104°F)
Storage temperatures	-40°C – + 70°C (-40°F – + 158°F)
Adjustable to head sizes in	30 increments
Color	Yellow
Environmental properties	Waterproof and UV-resistant
Head comfort	Foam pad
Tape liners	Easy pull tape liner
Headstrap	Medical grade adhesive



Ambu® Head Wedge®

MATERIAL


Product	Material
Head wedge base	Polypropylene
Actuator	Polyamid (PA) (Nylon)
Foam pad	Polyethylene (Recyclable PE)
Headstrap	Polyester (UP)




ORDERING SPECIFICATIONS

Item no.	Description
000264034	Cervical Immobilization Device w/ Head Strap, 10/case
000264036	Cervical Immobilization Device w/ Head 2 Strap, 10/case

Head Wedge can also be purchased with Perfit ACE and Mini Perfit ACE



Ambu® Head Wedge®



Ambu® Perfit ACE®
Ambu® Mini Perfit ACE®

Ambu FOREVER FORWARD

Ambu, Inc.
6721 Columbia Gateway Drive, Suite 200
Columbia, Maryland 21046
Tel. 800 262 8462
Fax 800 262 8673
ambuUSA.com

Questions? Contact your Ambu representative
at **800-262-8462** or visit **www.ambuUSA.com**

CE 2797 US: Rx only