

Ambu® Junior

DATASHEET

The Ambu Junior is a lifelike and cost-efficient training manikin that is designed to satisfy all training requirements.

The Ambu Junior has the appearance of a 3-year-old child with a body size of 34 inches.

The lifelike anatomy of the manikin helps to determine the correct hand placement when performing chest compressions.

With Ambu Junior, optimal training of chest compressions and ventilations are achieved. The Manikin offers realistic ventilations including ventilation resistance and expiration capability. The unique hygienic system prevents cross-contamination.

The closed-shaped torso, with its anatomically correct landmarks, allows students to learn and train on the use of Automatic External Defibrillators (AED) and the placement of pads.

Realistic Ventilation

The students will feel the realistic airway resistance and see the chest wall start rising when the head is correctly tilted with proper ventilations.

Just like in real life, Ambu Junior provides feedback when students look, listen and feel for ventilation.

Water Rescue

The Ambu Junior is waterproof, enabling CPR training for drowning victims.



KEY BENEFITS

Cost-efficient manikin for training of BLS algorithms

Lifelike anatomy of a child

Realistic training of chest compressions

Ambu hygienic system for mouth-to-mouth ventilation

Suitable for bag-mask ventilation

Visible chest rise

Airway only opens in correct head position

Realistic placement of AED defibrillation pads for training the use of AED

Designed for water rescue

Ambu

SPECIFICATIONS

Weight	4.1 kg
Dimension (length x width x height)	86.4 x 28.0 x 17.8 cm

ORDERING SPECIFICATIONS

Item No.	Description
A 282 010 000	Ambu® Junior incl. 5 mouth-nose pieces, 5 sealing rings, 25 head bags, T-Shirt, pants and carry bag

Accessories

100-2082U	Ambu Junior mouth-nose piece (set of 10)
282 010 501	Ambu Junior head bag (set of 100)
282 010 502	Ambu Junior sealing ring (set of 10)



Ambu Hygienic System



Ambu Junior Head

Questions? Contact your Ambu representative at
800-262-8462 or visit [www.ambuUSA.com/
ambu-junior](http://www.ambuUSA.com/ambu-junior)

Ambu FOREVER FORWARD

Ambu Inc.
6721 Columbia Gateway Drive
Suite 200
Columbia, MD 21046
Tel. 800 262 8462
Fax 800 262 8673
ambuUSA.com

CE