

BRONCHOSCOPY LIKE YOU'VE NEVER SEEN IT BEFORE

A NEW ERA FOR SINGLE-USE BRONCHOSCOPY



Ambu® aScope™ 5 Broncho

Ambu



ONE SOLUTION FOR A BROAD RANGE OF PROCEDURES

With aScope 5 Broncho, you're not limited by the number of available endoscopes. You have a variety of sterile bronchoscopes in different sizes available whenever you need them for diagnostic, therapeutic and peripheral bronchoscopy procedures.

Diagnostic and therapeutic procedures

The aScope 5 Broncho HD 5.0/2.2 and 5.6/2.8 endoscopes offer excellent maneuverability, even when working with endotherapy instruments, including high-frequency tools. This, combined with a bright, clear image, makes them suitable for a wide range of diagnostic and therapeutic procedures.



Ambu reinvents specimen collection

The single-use Ambu® aScope™ 5 Broncho Sampler Set is a unique integrated sampling solution that protects the sample, simplifies workflow and improves safety.

It includes 2x 60mL sample container to address the need for larger samples in BAL procedures.



The Ambu® aScope™ 5 Broncho Sampler Set

Protects the sample from start to finish

- Closed-loop system reduces the risk of sample contamination and loss common with Lukens traps
- Minimizes exposure points
- Protects the HCP from exposure to airway secretions

Simple to operate without assistance

- Seamlessly integrated with the aScope 5 Broncho endoscope
- Allows you to switch between suction and sampling at the turn of a knob, with no need to switch tubing
- High degree of control when navigating the lungs

Sterile and ready to go

- Single-use solution is sterile straight from the package
- Always available when you need it
- Clicks right onto the aScope 5 Broncho endoscope



Peripheral bronchoscopy procedures

aScope 5 Broncho 4.2/2.2 offers a unique opportunity to take advantage of a regular-sized working channel in a thin bronchoscope. The thin size and high bending angles of 210°/210° enable excellent maneuverability.

What's more, it is compatible with commonly used endotherapy instruments and active tools. These capabilities make aScope 5 Broncho 4.2/2.2 ideal for peripheral bronchoscopy procedures.

A single-use ultrathin option

With its thin outer diameter, aScope 5 Broncho 2.7/1.2 can be used in smaller airways without compromising maneuverability, image quality or access to the working channel. In fact, they offer a 120° field of view and 3-100 mm depth of field. The ultrathin aScope 5 Broncho 2.7/1.2 is suitable for ET tubes as small as 3.2.



ALL THE STERILITY AND EFFICIENCY OF SINGLE-USE

- **Sterile:** No risk of patient-to-patient contamination
- **Always available:** No reprocessing or repair bottlenecks
- **Consistent performance:** No wear and tear affecting the bending performance or imaging



PERFORMANCE COMPARABLE TO REUSABLES

- **High performance:** Excellent maneuverability and high-quality imaging
- **Greater flexibility:** A variety of readily available bronchoscopes in 4 different sizes
- **A unique sampling solution:** The aScope 5 Broncho Sampler Set simplifies BAL procedures (see other side)

Ambu® aScope™ 5 Broncho

BRONCHOSCOPY LIKE YOU'VE NEVER SEEN IT BEFORE

A new era for single-use bronchoscopy

aScope™ 5 Broncho flexible single-use bronchoscope

Item code	Product	Quantity
622001000US	Ambu® aScope™ 5 Broncho HD 5.6/2.8	5
621001000US	Ambu® aScope™ 5 Broncho HD 5.0/2.2	5
620001000US	Ambu® aScope™ 5 Broncho 4.2/2.2	5
624001000US	Ambu® aScope™ 5 Broncho 2.7/1.2	5

aScope™ 5 Broncho Sampler Set flexible single-use bronchoscope and specimen sampling system

Item code	Product	Quantity
622002000US	Ambu® aScope™ 5 Broncho HD 5.6/2.8 Sampler Set	5
621002000US	Ambu® aScope™ 5 Broncho HD 5.0/2.2 Sampler Set	5

Single-use sample container

Item code	Product	Quantity
480 500 000US	Ambu® aScope™ BronchoSampler 60 SC	10

Displaying and processing units

Item code	Product	Quantity
505001000US2	Ambu® aBox™ 2	1
405011000	Ambu® aView™ 2 Advance	1

aCart™

Item code	Product	Quantity
KU.9713.800	Ambu® aCart™ Compact KU	1
HA.3000.991	Ambu® aCart™ Compact Vario Mount HA	1



Ambu, Inc.
6721 Columbia Gateway Drive
Suite 200
Columbia, MD 21046
Tel. 800 262 8462
Fax 800 262 8673
ambuUSA.com

US Rx Only

PUB - 001763 - V01 - 08/2023 - Ambu USA. Technical data may be modified without further notice.