

# Ambu® aScope™ 4 RhinoLaryngo Slim

## DATASHEET

### SINGLE-USE FLEXIBLE RHINOLARYNGOSCOPE FOR FEES EXAMS

#### SHARP IMAGE DELIVERY

Provides a clear view of anatomical structures and simplifies navigation of the upper airway

#### ERGONOMIC, LIGHTWEIGHT DESIGN

Maximizes comfort when performing lengthier swallow assessments

#### SLIM SCOPE DIAMETER

Designed for increased patient comfort during exams

#### ALWAYS-READY, SINGLE-USE PLATFORM

Optimizes scheduling flexibility and improves workflow with no cleaning or repair necessary



### KEY BENEFITS FOR SPEECH-PATHOLOGY

Provides immediate availability - ready when needed to perform swallow study FEES exams


Eliminates reprocessing after exams, along with point-of-care difficulties and endoscope transport.

Cuts out capital or service/repair expenses after repeated usage of traditional endoscopes

Transition from fiber-optic/camera setups and bulky, stationary video towers

Connects to highly portable, compact HD display that is plug and play, mounts easily to IV poles, and operates on battery power when needed

Nation-wide network of Clinical Training Specialists to support new users in training and hospital/clinic installation



130°/130° bending angle allows for easy viewing of hard-to-access anatomical structures. When combined with the portable 12.8" HD aView 2 Advance monitor, this single-use platform provides an effective, easy-to-use solution for swallow assessments.

# Ambu

## SPECIFICATIONS

### Optical system

Field of view	85°
Direction of view	0° (forward viewing)
Depth of field	6 - 50 mm
Illumination method	LED

### Insertion portion

Bending section	130° up, 130° down
Insertion cord diameter	3.0 mm (0.12")
Distal tip diameter	3.5 mm (0.14")
Maximum diameter of insertion portion	3.5 mm (0.14")
Working length	300 mm (11.8")

### Operating environment, storage & transportation

Temperature - transportation	10 - 40° C (50 to 104° F)
Temperature - operating environment	10 - 40° C (50 to 104° F)
Temperature - storage	10 - 25° C (50 to 77° F)
Relative humidity	30 - 85%
Atmospheric pressure	80 - 109 kPa

### Sterilization


Method of sterilization	ETO
-------------------------	-----

## ORDERING SPECIFICATIONS\*


Item no.	Quantity
----------	----------

510001000 - aScope 4 RhinoLaryngo Slim	5
--	---

\* aScope 4 RhinoLaryngo Slim is not available for sale in all countries. Please contact your local Ambu sales representative.



**Ambu® aScope™ 4 RhinoLaryngo Slim  
and Ambu aView™ 2 Advance  
on Ambu aCart™ Compact**



**Ambu** FOREVER FORWARD

Ambu, Inc.  
6721 Columbia Gateway Drive, Suite 200  
Columbia, Maryland 21046  
Tel. 800 262 8462  
Fax 800 262 8673  
ambuUSA.com

**CE** 2797 US: Rx only

LBL - 008714 - V03 - 04/2025 - Ambu USA. Technical data may be modified without further notice.