

aCart™ Compact

DATASHEET

Designed for convenience

For teams who perform endoscopic procedures, aCart Compact is a modern, mobile workstation that adds value to the Ambu single-use product portfolio. The aCart Compact workstation provides convenient storage and portability without taking up a lot of space.

A modern mobile workstation

The aCart Compact workstation enables you to quickly and seamlessly respond to an emergency or change clinical environments. It is also ideal for clinical use where space is limited.

Robust and portable

The design and choice of materials make the aCart Compact workstation light enough to easily transport, yet robust and stable enough to protect your valuable endoscopic equipment.

Everything is within reach

The aCart Compact workstation puts everything you need within reach thanks to a:

- Storage hook for spare single-use endoscopes
- Basket for accessories, cables or airway devices
- Preparation shelf
- Dual scope holder for single-use endoscopes that have been prepared for use and that are in use
- VESA monitor mount, which the Ambu displaying and processing units can be attached to
- Adjustable IV pole
- Cable duct for efficient cable management



TECHNICAL SPECIFICATIONS

Dimensions

Width	21 in.
Height without IV pole	52 in.
Height with IV pole min./max.	60/93 in.
Depth	24 in.
Package Dimensions W/D/H	46/23/9 in.

Weight

Total	51 lbs.
-------	---------

Load Capacity

Maximum load	77 lbs.
Shelve	22 lbs.
Monitor holder	31 lbs.
Basket	7 lbs.
IV Pole	11 lbs.
Optional Vario Mount	18 lbs.

Device Classification

Europe	Class I, CE certified according to DIN EN 60601-1
USA	Class I, 510(k) exempt from premarket notification

Operating environment

Operating temperature	50 - 104°F
Operating relative humidity	30 - 75%
Operating atmospheric pressure	70 - 106kPa

Storage environment

Storage temperature	-13 - 158°F
Storage relative humidity	10 - 95%
Storage atmospheric pressure	50 - 120kPa

ORDERING SPECIFICATIONS

Description	Reference number
aCart Compact	KU.9713.800
aCart Compact Vario Mount	HA.3000.991

Have questions about the aCart™ Compact workstation?

Please contact your local Ambu representative at **800-262-8462** or visit **www.ambu.USA.com**

Ambu FOREVER FORWARD

Ambu, Inc.
6721 Columbia Gateway Drive, Suite 200
Columbia, MD 21046
Tel. 800 262 8462
Fax 800 262 8673
ambuUSA.com