

aCart™ Plus

DATASHEET

Designed for stability and convenience

For teams who perform procedures in the endoscopy suite, the aCart Plus workstation is designed for storage and convenience.

Robust and efficient

The robust design makes the aCart Plus workstation durable enough to meet the demands of everyday use. The materials and surface finishes are high quality.

Everything is within reach

The smart layout of the cart puts the tools you need within reach with:

- Ample shelf space for medical devices
- Medical rail for attachable accessories
- Sterile water bottle and CO₂ cylinder holders
- VESA mount for up to 32" monitor
- Scope holder for single-use endoscopes that have been prepared for use and that are in use
- Two cable ducts for efficient cable management



TECHNICAL SPECIFICATIONS

Dimensions

Width	38 in. (total); 24 in. (without CO2 cylinder holders)
Height without scope holder	68 in.
Height w/ scope holder min./max.	78/112 in.
Depth	25 in.
Package Dimensions W/D/H	28/28/68 in.

Weight

Total	95 lbs.
-------	---------

Load Capacity

Maximum load	331 lbs.
Shelve (each)	110 lbs.
Monitor holder	31 lbs.
CO ₂ cylinder holder (each)	44 lbs.
Endoscope holder	11 lbs.
Medical Accessory Rail	22 lbs.

Device Classification

Europe	Class I, CE certified according to DIN EN 60601-1
USA	Class I, 510(k) exempt from premarket notification

Operating environment

Operating temperature	50 - 104°F
Operating relative humidity	30 - 75%
Operating atmospheric pressure	70 - 106kPa

Storage environment

Storage temperature	-13 - 158°F
Storage relative humidity	10 - 95%
Storage atmospheric pressure	50 - 120kPa

ORDERING SPECIFICATIONS

Description	Reference number
aCart Plus	KU.9740.800

Have questions about the aCart™ Plus workstation?

Please contact your local Ambu representative at **800-262-8462** or visit **www.ambu.USA.com**



Ambu, Inc.
6721 Columbia Gateway Drive, Suite 200
Columbia, MD 21046
Tel. 800 262 8462
Fax 800 262 8673
ambuUSA.com