



Ives MR Conditional EEG Electrode System

Head Diagram
 Disposable Black Conductive
 Plastic Electrode Set
 Model: DCPE-21/B
 (Ambu-21)

Conductive Plastic
 Electrode: Black Cup
 Version w Silver, no label



MR Conditional
 See user manual 



LEFT (GRN)
 Fp1 BRN
 F3 RED
 C3 ORG
 P3 YEL
 O1 GRN
 F7 BLU
 T3 PUR
 T5 GRY
 Cz WHT
 GND BLK

RIGHT (RED)
 Fp2 BRN
 F4 RED
 C4 ORG
 P4 YEL
 O2 GRN
 F8 BLU
 T4 PUR
 T6 GRY
 REF WHT
 Fz BLK
 Pz PNK

J.R. Ives
 December 22, 2017