



Ives MR Conditional EEG Electrode System

Head Diagram
 Disposable
 Conductive Plastic Electrode
 Set
 Model: DCPE-25/B
 (Ambu-25)



MR Conditional
 See user manual



LEFT (GRN)
 Fp1 BRN
 F3 RED
 C3 ORG
 P3 YEL
 O1 GRN
 F7 BLU
 T3 PUR
 T5 GRY
 T2 WHT
 A1 BLK

RIGHT (RED)
 Fp2 BRN
 F4 RED
 C4 ORG
 P4 YEL
 O2 GRN
 F8 BLU
 T4 PUR
 T6 GRY
 T1 WHT
 A2 BLK

CENTRAL (YEL)
 Fz BRN
 Cz RED
 Pz ORG
 GND BLK
 REF WHT