

The Ambu® RDIC Resuscitator

Ambu is pleased to announce the availability of the world's only effective, truly portable manual resuscitation system designed and tested specifically for use in chemical or biological conditions



Main Features Include:

- Developed according to US military specifications
- Butyl outer cover for optimal resistance to penetration during ventilation in toxic atmospheres
- Can be used in all conventional ventilation emergencies
- NATIONAL STOCK NUMBER (NSN): #6515-01-338-6602
- PART NUMBER: #233-005-000

THE AMBU RDIC IS SUPPLIED AS A COMPLETE SYSTEM INCLUDING:

- A genuine ***Ambu Mark III*** reusable resuscitator which is equipped with impermeable elastomeric components on the compression bag and patient valve (butyl rubber)
- A highly flexible butyl rubber extension tube which offers unparalleled flexibility during casualty transport
- A contoured face mask with clear dome which permits better monitoring of the patient
- A specialized inlet filtration system, including C2 canister, for filtering contaminated air
- A durable, high-impact molded carry case

Ambu 

Ambu Inc. 6740 Baymeadow Drive, Glen Burnie, MD 21060
Phone: (800) 262-8462 Fax: (800) 262-8673
www.ambuUSA.com