

ONE SOLUTION FOR ALL YOUR AIRWAY VISUALIZATION NEEDS



Ambu

Ambu® SureSight® Video Laryngoscopy System

ONE SOLUTION FOR ALL YOUR AIRWAY VISUALIZATION NEEDS

The Ambu airway visualization platform is designed to save space and reduce task complexity. It is an integrated platform that puts your key airway visualization tools within reach and ready to use.

The Ambu airway visualization platform is the only complete airway visualization platform featuring a one-lung ventilation solution, including all of the devices recommended for intubation and airway assessment*.

*Recommended by the UK Difficult Airway Society; American Society for Anesthesiologists; Australian and New Zealand College of Anaesthetists & Faculty of Pain Medicine



NEW!

Ambu® SureSight® Video Laryngoscope

Intubation - Maximize your focus on the patient



Ambu® aScope™ 5 Broncho

Bronchoscopy - High performance and flexibility

Ambu® aScope™ 5 Broncho Sampler Set

A unique integrated BAL/BW solution



Ambu® aScope™ 4 Broncho

Bronchoscopy - Portability and essential performance

Ambu® aScope™ 4 Broncho Sampler Set

A unique integrated BAL/BW solution



Ambu® aView™ 2 Advance Ambu® aBox™ 2

Endoscopy systems - Your control center for key devices



Ambu® VivaSight™ 2 DLT

OLV procedures - Built-in camera increases patient safety



Ambu® aScope™ 4 RhinoLaryngo

Rhinolaryngoscopy - For full examination of the upper airway in difficult airway cases



Together with the visualization devices, Ambu laryngeal masks, face masks and resuscitators, Ambu comprises a comprehensive airway management solution.

The Ambu® SureSight® Connect Video Laryngoscope Solution

Ambu® aCart™ Compact


Everything within reach and ready to go. A compact, mobile workstation for the SureSight Connect video laryngoscope solution and your other key devices



Comprehensive blade portfolio

VL blades for both routine and difficult intubations as well as emergent or planned procedures.

- Miller 0
- Miller 1
- MAC 3
- MAC 4
- Hyperangulated 1
- Hyperangulated 2
- Hyperangulated 2 Channeled
- Hyperangulated 3
- Hyperangulated 3 Channeled
- Hyperangulated 4



Ambu® aView™ 2 Advance Ambu® aBox™ 2 (not shown here)

Plug-and-play system that serves as the control center for key Ambu airway visualization devices

Ambu® SureSight® Connect

An intuitive, ergonomic and reliable video laryngoscope

- Audible click ensures that the blade is properly attached
- Ergonomic design for secure grip and optimal control
- Anti-fog lens for a clear and consistent image
- Easy-to-reach video/camera button

EVERYTHING YOU NEED WITHIN REACH AND READY TO GO

The Ambu® SureSight® Connect video laryngoscope is part of our trusted airway visualization platform. Having one uniform, intuitive, and familiar system helps you stay focused on the patient as you switch between devices and move through the airway management algorithm.

One system, multiple visualization devices

With the Dual View feature in the Ambu endoscopy system, you are able to simultaneously use two of our airway visualization devices without needing separate screens for difficult airways.

Streamlines work with patient imaging data

The Ambu endoscopy system offers connectivity with Electronic Health Records (EHR) to simplify the important task of sharing and transferring patient data.

[Learn more about the Ambu intubation solution.](#)

THE WORLD'S MOST TRUSTED SINGLE-USE BRONCHOSCOPY SOLUTION

It has been tried and tested in millions of cases, at thousands of hospitals around the world. It also provides you with a unique, integrated BAL/BW sampling set.

aScope 4 Broncho: Portable, lightweight, easy to set up

Offers the portability and essential performance needed for a variety of procedures in the OR to routine bronchoscopy in the ICU.

aScope 5 Broncho: A new standard for single-use bronchoscopes

Combines the high level of maneuverability and imaging needed for advanced diagnostic and therapeutic procedures with the sterility and efficiency of the single-use concept.

aScope Sampler Set: Designed for superior lavage

Reduces the risk of sample contamination and loss common with Lukens traps, protecting the sample from start to finish.

[Learn more about the Ambu bronchoscopy solution.](#)

The Ambu® aScope™ Broncho Solution

The Ambu® aScope™ Broncho portfolio provides two single-use endoscopy solutions with a variety of sizes available in storage whenever you need them.

Ambu® aView™ 2 Advance Ambu® aBox™ 2 (not shown here)

Plug-and-play system that serves as the control center for key Ambu airway visualization devices

Ambu® aCart™ Compact

Everything within reach and ready to go. A compact, mobile workstation for the SureSight Connect video laryngoscope solution and your other key devices

Ambu® aScope™ 5 Broncho

- 2.7/1.2 An ultrathin bronchoscope suitable for ET tubes as small as 3.5 mm
- 4.2/2.2 A regular-sized working channel in a thin bronchoscope
- 5.0/2.2 and 5.6/2.8 Diagnostic and therapeutic scopes for advanced bronchoscopies



Ambu® aScope™ 4 Broncho

- 3.8/1.2 From difficult intubation to placement and control of DLTs and bronchial blockers in the OR
- 5.0/2.2 From intubation and training to bedside bronchoscopy, including PDT procedures
- 5.8/2.2 With a suctioning channel of 2.8 mm, it can be used for secretion management



Ambu® aScope™ Broncho Sampler Sets *A closed-loop single-use system*

- Ambu® aScope™ 4 Broncho Sampler Set
- Ambu® aScope™ 5 Broncho Sampler Set

The Ambu® VivaSight™ 2 OLV Solution

Ambu® aCart™ Compact

Everything within reach and ready to go. A compact, mobile workstation for the SureSight Connect video laryngoscope solution and your other key devices

Ambu® VivaSight™ 2 DLT

Built-in camera at the tip of the tracheal lumen enables continuous visualization

Ambu® aView™ 2 Advance Gen. 2

Plug-and-play system that serves as the control center for key Ambu airway visualization devices

Ambu® aScope™ 4 and 5 Broncho

A variety of different sizes are always available to handle the various situations that can arise during OLV procedures



CONTINUOUS VISUALIZATION INCREASES PATIENT SAFETY AND IMPROVES WORKFLOW

The Ambu® VivaSight™ 2 portfolio offers a safe and cost-effective way to perform One-Lung Ventilation (OLV). Although similar to conventional DLT and SLT solutions, VivaSight 2 has a built-in camera at the tip of the tracheal lumen. This enables you to immediately detect and correct tube malpositioning and dislocation without compromising ventilation of the patient.

No need to set up a second system

VivaSight 2 significantly reduces the need to use a bronchoscope. In cases where you do need a bronchoscope, however, you can take advantage of Dual View with Ambu® aScope™ 4 or 5 Broncho and aView 2 Advance (Gen. 2). There's no need to set up a second system.

[Learn more about the Ambu VivaSight 2 OLV solution.](#)

HELPS YOU STREAMLINE YOUR WORKFLOW AND PERFORM PROCEDURES CONFIDENTLY

The Ambu® aScope™ 4 RhinoLaryngo solution is a single-use endoscopy solution that offers excellent imaging, is always available and easy to take with you.

Excellent imaging with integrated video and audio

- Easily identify anatomical structures with high endoscope image quality and the Ambu full-HD endoscopy system
- Dictate notes as you perform the procedure with synchronized audio and video
- Review recordings of fast movements more carefully with slow-motion video playback
- Isolate and save key frames for reporting, training and to get second opinions

[Learn more about the Ambu RhinoLaryngoscope solution.](#)

The Ambu® aScope™ RhinoLaryngo Solution

The Ambu® aScope™ RhinoLaryngo portfolio provides two single-use endoscopy solutions available in storage whenever you need them.

Ambu® aView™ 2 Advance

Plug-and-play system that serves as the control center for key Ambu airway visualization devices




Ambu® aScope™ 4 RhinoLaryngo

- 3.0 outer diameter: An excellent choice for nasal endoscopies, laryngoscopies, and swallow assessments
- 5.0/2.0 An excellent choice for therapeutic procedures



Ambu® aCart™ Compact

Everything within reach and ready to go. A compact, mobile workstation for the SureSight Connect video laryngoscope solution and your other key devices



THE AMBU AIRWAY VISUALIZATION PLATFORM

One single-use platform; multiple long-term benefits

The Ambu endoscopy systems are used across many clinical areas. This means you can use the familiar Ambu endoscopy system already present in the department when you are called to establish an airway. Having one high-quality platform designed to work together means:

- A familiar and intuitive system
- One trustworthy supply partner to simplify procurement and provide ongoing support
- Access to a portfolio of single-use endoscopes and OLV devices, which can reduce costs associated with capital investments, reprocessing, and repairs

ORDERING SPECIFICATIONS

Ambu® SureSight® Video Laryngoscope

Item code	Product	Quantity
390511000	Ambu SureSight Connect US	1
390101000	Ambu SureSight Blade Miller Size 0 US	5
390111000	Ambu SureSight Blade Miller Size 1 US	5
390211000	Ambu SureSight Blade Hyperangulated Size 1 US	5
390221000	Ambu SureSight Blade Hyperangulated Size 2 US	5
390521000	Ambu SureSight Blade Hyperangulated Size 2 Channeled US	5
390231000	Ambu SureSight Blade Hyperangulated Size 3 US	10
390531000	Ambu SureSight Blade Hyperangulated Size 3 Channeled US	10
390241000	Ambu SureSight Blade Hyperangulated Size 4 US	10
390331000	Ambu SureSight Blade MAC Size 3 US	10
390341000	Ambu SureSight Blade MAC Size 4 US	10

Ambu® aScope™ 5 Broncho / Ambu® aScope 5 Broncho™ Sampler Set

Item code	Product	Quantity
621001000US	Ambu® aScope™ 5 Broncho HD 5.0/2.2	5
622001000US	Ambu® aScope™ 5 Broncho HD 5.6/2.8	5
620001000US	Ambu® aScope™ 5 Broncho 4.2/2.2	5
624001000US	Ambu® aScope™ 5 Broncho 2.7/1.2	5
621002000US	Ambu® aScope™ 5 Broncho HD 5.0/2.2 Sampler Set	5
622002000US	Ambu® aScope™ 5 Broncho HD 5.6/2.8 Sampler Set	5
620002000US	Ambu® aScope™ 5 Broncho 4.2/2.2 Sampler Set	5

Ambu® aScope™ 4 Broncho / Ambu® aScope 4 Broncho™ Sampler Set

Item code	Product	Quantity
476101000	Ambu® aScope™ 4 Broncho Slim	5
477101000	Ambu® aScope™ 4 Broncho Regular	5
478101000	Ambu® aScope™ 4 Broncho Large	5
477501000	Ambu® aScope™ 4 Broncho Regular Sampler Set	5
478501000	Ambu® aScope™ 4 Broncho Large Sampler Set	5

Ambu® VivaSight™ 2

Item code	Product	Quantity
412351000	Ambu® VivaSight™ 2 DLT Kit 35 Fr - LEFT	1
412371000	Ambu® VivaSight™ 2 DLT Kit 37 Fr - LEFT	1
412391000	Ambu® VivaSight™ 2 DLT Kit 39 Fr - LEFT	1
412411000	Ambu® VivaSight™ 2 DLT Kit 41 Fr - LEFT	1

Ambu® Display Units and Carts

Item code	Product	Quantity
405011000US2	Ambu® aView™ 2 Advance Gen. 2 Displaying and Processing Unit	1
505001000	Ambu® aBox™ 2 Displaying and Processing Unit	1
KU.9713.800	Ambu® aCart™ Compact	1
KN.1215.991	Ambu® aCart™ Compact Vario Mount	1
KN.1245.800	Ambu® aCart™ Compact Universal Endoscope Holder	1

Ambu® aScope™ 4 RhinoLaryngo

Item code	Product	Quantity
510101000	Ambu® aScope™ 4 RhinoLaryngo Slim	5
512101000	Ambu® aScope™ 4 RhinoLaryngo Intervention	5

SureSight Photo Gallery

THE AMBU COMMITMENT

Ambu is on a mission to deliver airway management solutions that help you care for patients and save lives. At the same time, as a single-use medical device company, we are acutely aware of our obligation to work toward more sustainable solutions.

Taking leaps forward in sustainability

Advancing in the area of sustainability is one of the four strategic pillars of our company. This strategy focuses on circular design principles that help us lower our products' carbon footprint and activities throughout our value chain that help us work toward our net zero emissions goal.

One concrete example of this is our pioneering use of bioplastic material, which is a mix of fossil-based and 2nd generation bio feedstock, in our endoscope handles, laryngeal mask cuff protector packaging and other products.

By year end 2024, all our endoscopes were made with bioplastic material in the handles. In the future, we plan to introduce the material in other parts of our endoscopes. In addition, SureSight is phthalate and PVC-free.

Contact us for more information

Learn more about the Ambu airway visualization platform, including solutions for intubation, bronchoscopy and One Lung Ventilation.

Questions? Contact your Ambu representative at **800-262-8462** or visit AmbuUSA.com/SureSight

Ambu, Inc.
6721 Columbia Gateway Drive
Suite 200
Columbia, MD 21046
Tel. 800 262 8462
Fax 800 262 8673
AmbuUSA.com

US: Rx Only

PUB-002088 - V02 - 05/2025 - Ambu USA - Technical data may be modified without further notice.