

Ambu[®] Res-Cue Key



Key Benefits

- Semi-transparent face shield and one-way valve allowing rescuer to visually check for patient's lip color and vomitus
- Prevents contact with patient's saliva or blood
- Suitable for both adults and children
- Elastic ear bands fit over patient's ears to keep the mask in place
- Can be carried on a key chain, in a pocket, a small handbag and a first-aid kit

Ambu[®] Res-Cue Key

The Ambu Res-Cue Key is a proven, easy to use, and inexpensive barrier with a one-way valve that prevents direct contact with patient's mouth, nose and face.

The Ambu Res-Cue Key helps protect rescuer and increases bystanders' motivation for CPR in emergency situations.

Standard resuscitation techniques can be used without modification.

Comes with Ambu logo or as private labelling. Available in hardcase, softcase or plastic bag.

Ordering specifications

Ambu Res-Cue Key	
248 201 102	Res-Cue key in Softcase with Ambu logo Black
248 201 103	Res-Cue key in Softcase with Ambu logo Blue
248 201 104	Res-Cue key in Softcase with Ambu logo Red
Ambu Res-Cue Key private label	
248 201 303	Res-Cue key in Softcase with private label Blue
248 201 304	Res-Cue key in Softcase with private label Red
248 201 302	Res-Cue key in Softcase with private label Black

Materials


Product	Material
Mask	Polyethylene
Softcase	Nylon
Hardcase	Styrolux

Size


Type	
Folded	3 x 2 x 1 cm
Hardcase	3.7 x 3.7 x 1.5 cm
Softcase	5 x 5 x 1.5 cm
Weight app	10 g

Softcase

Printing area (logo size)	14 x 39 cm Available in colors: black, blue, red, yellow
Order	
Minimum	250 pieces



Ambu Res-Cue Key in use



Ambu Res-Cue Key private label