

Ambu® Single-Use Urology Portfolio

ONE UROLOGY SOLUTION FOR EVERY PROCEDURE

Ambu



ONE PLATFORM, ONE UROLOGY SOLUTION



aScope™ 4 Cysto


Developed by the world leader in single-use endoscopy, the aScope 4 Cysto is the first cystoscope to the Ambu urology portfolio. Ideal for routine cystoscopies, providing dependable performance with every procedure and minimal setup time.

- Boost physician efficiency and productivity
- Maximize operational efficiencies
- Allows for reallocation of valuable resources
- Enhanced ergonomics with easy setup and use
- Improved patient satisfaction

aScope™ 5 Cysto HD

Advanced HD precision for enhanced visualization, perfect for complex cases requiring the highest level of detail – all the features and capabilities of the aScope 4 Cysto, plus more.

- Advanced Red Contrast (ARC) technology
- Enhanced compatibility with electro-frequency tools
- Easy visualization of vascularity and suspicious lesions on the full HD display unit



aScope™ 5 Uretero

Versatile and efficient, offering reliable ureteroscopy performance with the flexibility to handle various patient needs.

- 7.9 Fr tapered distal tip designed to access difficult-to-reach areas
- No downtime and costs associated with the cleaning and repair of reusable ureteroscopes
- Lightweight, ergonomic design enables excellent torque and enhanced maneuverability

ONE PLATFORM, MAXIMUM VALUE

Ambu's endoscopy systems offer an intuitive solution designed for efficiency. Built on a common platform, they streamline workflow and provide an effective way to standardize across clinical areas.



aView 2 Advance: The portable solution provides a 12.8-inch, HD battery-powered display unit that fits into small spaces and can be mounted on an IV pole.

WorkList Software Updates

- Improved Workflow: Import and export patient data effortlessly with the new worklist feature, fully integrated with your EHR system.
- Enhanced Cybersecurity: Benefit from industry-leading WPA2 and WPA3 compliance for the most secure Wi-Fi connectivity.
- All-in-One Solution: Eliminate the need for third-party devices like smartphones or USBs – everything you need is built in.

aCart™ Compact is a modern, mobile workstation providing easy convenient storage for endoscopy procedures in a small footprint – convenient, portable, versatile and hygienic.



aBox 2: The OR solution offers a total endourology experience, powering both aScope 5 Uretero and aScope 5 Cysto HD with dual imaging capabilities all while maintaining a small integrated footprint.



LEADING THE INDUSTRY IN ECO-CONSCIOUS DESIGN

All of our endoscope handles are now made with bioplastics, which reduces the carbon footprint of our single-use endoscopes. This is just part of our push on several fronts to lead the way in developing sustainable practices that also help our customers reduce their environmental footprints:

- Use of bioplastic material reduces the carbon footprint of the ABS plastics used in the production of our ureteroscope handles by 70%
- Bioplastics innovation in device design and packaging cuts carbon emissions and enhances recyclability
- Our circular business model enables comprehensive product-level environmental reporting



Learn about Ambu's sustainability initiatives, including our innovative use of bioplastics.



ORDERING SPECIFICATIONS

Item No.	Description
601101000	aScope™ 4 Cysto, Standard Deflection
601001000 (Canada)	aScope™ 4 Cysto, Standard Deflection
600101000	aScope™ 4 Cysto, Reverse Deflection
600001000 (Canada)	aScope™ 4 Cysto, Reverse Deflection
603101000	aScope™ 5 Cysto HD, Standard Deflection
602101000	aScope™ 5 Cysto HD, Reverse Deflection
604101000	aScope™ 5 Uretero, Standard Deflection
605101000	aScope™ 5 Uretero, Reverse Deflection
505001000US2	aBox™ 2
405011000US2	aView™ 2 Advance (Gen. 2)
FLO 30	Hemostasis Valve FLO 30
K12-11412	Stopcock 4-Way. 50 psi (3.5 bar) Swivel Nut, Large Bore
K12-11413	Hemostasis Valve. AccessPLUS™, Large Bore
KU.9713.800	aCart™ Compact

Ambu FOREVER FORWARD

Ambu, Inc.
6721 Columbia Gateway Drive
Suite 200
Columbia, MD 21046
ambuUSA.com