

Ambu® aScope™ 4 Cysto

The single-use flexible digital cystoscope



Maximize Physician Efficiency

Scheduling and workflow flexibility to increase daily patient procedures – including last-minute add-ons

Maximize Staff Productivity

Reallocate staff resources from time-consuming reprocessing to patient-focused care

Budget Optimization

No repairs or reprocessing costs, and limited capital expenditures while increasing procedural capacity

Simplified Compliance

Sterile, patient-ready cystoscope without the complexities of reprocessing and monitoring

Patient Experience

Expand procedure locations to additional rooms or satellite clinics, and enable image/video review with patients



Developed by the world leader in single-use endoscopy, the aScope 4 Cysto is the newest addition to Ambu's line of endoscopes. When combined with the portable 12.8" HD aView2 Advance monitor, this single-use platform provides an effective solution for flexible cystoscopy procedures across various clinical settings.

Ambu® aScope™ 4 Cysto | Ready When Needed

Specifications

Optical system	
Field of view	120°
Depth of field	3 - 100 mm
Illumination method	LED
Insertion portion	
Bending angle*	210° up, 120° down
Insertion cord diameter	16.2 Fr / 5.4 mm (0.21")
Distal end diameter	16.2 Fr / 5.4 mm (0.21")
Maximum diameter of insertion portion	18 Fr / 6.0 mm (0.24")
Working length	390 mm (15.4")
Working channel	
Minimum instrument channel width**	6.6 Fr / 2.2 mm (0.086")
Operating environment, storage & transportation	
Temperature - transportation	10 - 40 °C (50 to 104 °F)
Temperature - operating environment	10 - 40 °C (50 to 104 °F)
Temperature - storage***	10 - 25 °C (50 to 77 °F)
Relative humidity	30 - 85%
Atmospheric pressure	80 - 109 kPa
Sterilization	
Method of sterilization	ETO
Material	
DEHP phthalates	Does not contain any DEHP phthalates
No allergenic substances	Not made with natural rubber latex
Dimensions	
Weight	158 g






aScope 4 Cysto offers consistent quality every time



High bending angles of 210°/120°

Order specifications

Product	Item no.	Quantity
Ambu® aScope™ 4 Cysto, Standard Deflection	601001000	5 pcs
Ambu® aScope™ 4 Cysto, Reverse Deflection	600001000	5 pcs
Accessories		
Stopcock 4-Way. 50 psi (3.5 bar). Swivel Nut, Large Bore****	K12-11412	50 pcs
Hemostasis Valve. AccessPLUS™, Large Bore****	K12-11413	25 pcs



Ambu A/S

Ambu Inc.
6230 Old Dobbin Lane
Columbia, MD 21045
Tel. 800 262 8462
Fax 800 262 8673
www.ambuUSA.com

* Please be aware that the bending angle can be affected if the insertion cord is not kept straight or have inserted endoscopic accessories.
 ** There is no guarantee that endoscopic accessories selected solely using this minimum working channel width will be compatible in combination.
 *** Storage under higher temperatures may impact shelf life.
 **** Not available in all countries.

US: Rx only

