

Ambu® aScope™ Duodeno

Sterile. Single-use. Problem solved.



Key Benefits

- Sterile
- No risk of patient cross-contamination
- Cost effective
- Eliminates the risk and cost of reprocessing
- Familiar design and function to that of conventional duodenoscopes





Distal tip showing elevator in use with a sphincterotome



Elevator control lever to raise (up) and lower (down) the elevator

Sterile straight from the pack, Ambu aScope Duodeno is a single-use endoscope with side viewing optics, deflectable tip and an elevator to control the position of compatible endoscopic accessories.

With aScope™ Duodeno, you have a new scope with reliable performance every time, and there is no risk of patient cross-contamination. What's more, there's no need for costly reprocessing or repair.



Specifications - Ambu® aScope™ Duodeno

Optical System	
Direction of View	Retro Side Viewing 6°
Field of View	130°
Insertion Section	
Distal End Outer Diameter	φ13.7 mm
Insertion Tube Outer Diameter	φ11.3 mm
Working Length	1240 mm
Instrument Working Channel	
Inner Diameter	φ4.2 mm (12.6 FR)
Bending Section	
Angulation Range	Up 120° / Down 90° Left 90° / Right 110°
Illumination Method	Dual LEDs
Dimensions	
Weight	700 g

Ordering Specifications

aScope Duodeno	Quantity
482001000	Packs of 6
Used with Ambu aBox Duodeno - 485001000US	

Ambu A/S

Baltorpbakken 13
 DK - 2750 Ballerup
 Denmark
 Tel. +45 7225 2000
 Fax +45 7225 2053
 www.ambu.com