

KingJackson Rees & KingMapleson

non-rebreathing circuit kits



Flexible. Control. Compatible.

King Systems' disposable modified Mapleson and Jackson-Rees non-rebreathing circuits offer the widest variety of designs and options available to provide convenient solutions for difficult non-rebreathing circuit configurations

Flexible.

- Multiple port and elbow options allow practitioners to customize their configuration.
- Standard and ultra-flex tubing offer practitioners complete flexibility during procedures.

Control.

- Adjustable thumb slide and tail valves offer total control over the delivery of medical gases.

Compatible.

- Compatibility with King Systems' breathing bags and face masks create a complete system.

KingJackson Rees & KingMapleson

non-rebreathing circuit kits

	Part Number	Case Quantity	10" 22mm Tubing	7' Oxygen Tubing	72" Fresh Gas Tubing	120" Fresh Gas Tubing	72" Fresh Gas Tubing	54" Bubble Line	96" Supply Line	Norman Elbow	Sample Port Elbow	Dual Port Elbow	Elbow	Hose Clip	Thumb Slide Valve	Tail Valve	Washington "R"	Standard Adaptor	*Breathing Bag	48" Scavenging Tubing and 19mm ID Connector
Jackson Rees																				
Non-Rebreathing Circuit Kit	60740	20		•						•					•			•		
Modification Non-Rebreathing Circuit Kit	60743Z	20		•						•					•					3L
Non-Rebreathing Circuit Kit	3701P	20	•		•					•										1 L Purple
Modification Non-Rebreathing Circuit Kit	3741P	20	•		•					•				•	•					1 L Purple •
Modification Non-Rebreathing Circuit Kit	3742B	20	•		•					•				•	•					2 L Blue •
Modification Non-Rebreathing Circuit Kit	3745Y	20	•		•					•				•	•					1/2 L Yellow •
Non-Rebreathing Circuit Kit	3771P-WSN	20	•		•					•										1 L Purple
Non-Rebreathing Circuit Kit	3771Y-WSN	20	•		•					•										1 L Yellow
Non-Rebreathing Circuit Kit	3772-WSNZ	20	•				•			•										2 L
Non-Rebreathing Circuit Kit	3773-WSNZ	20	•		•					•										3 L
Non-Rebreathing Circuit Kit	3801R-61L	20	•			•				•										1 L Red
Non-Rebreathing Circuit Kit	3841P-61	20	•		•					•				•	•					1 L Purple •
Non-Rebreathing Circuit Kit	3842B-61	20	•		•					•				•	•					2 L Blue •
Modification Non-Rebreathing Circuit Kit	3845Y-61	20	•		•					•				•	•					1/2 L Yellow •

*Breathing Bags are not made of Natural Rubber latex.

Mapleson C																				
Non-Rebreathing Circuit Kit	30840L	20							•				•		•					
Non-Rebreathing Circuit Kit	4211P	20		•							•				•					1 L Purple
Non-Rebreathing Circuit Kit	4211R	20		•							•				•					1 L Red
Non-Rebreathing Circuit Kit	4221B	20		•							•				•					2 L Purple
Non-Rebreathing Circuit Kit	4251Y	20		•							•				•					1/2 L Yellow

Sales Company
Ambu Inc.
 6740 Baymeadow Drive
 Glen Burnie, MD 21060
 Tel.: 1 410 768 6464
 1 800 262 8462
 Fax. 1 800 262 8673
 www.ambuUSA.com

Manufacturer
King Systems
 15011 Herriman Boulevard
 Noblesville, IN 46060
 Tel. 317 776 6823
 Tel. 800 642 5464
 Fax 317 776 6827
 www.Kingsystems.com

US: Rx only