

**Ambu**

**THE SINGLE-USE  
SOLUTION FOR  
CYSTOSCOPY**

Reimagine Procedure Days

---

aScope™ 4 Cysto





**A flexible video cystoscopy system for office-based urology practices to improve workflow, staff efficiency, and patient throughput without the need for significant capital investment, reprocessing, or repairs.**

Eliminate traditional expenses and equipment challenges in your practice

- No bulky towers
- No light sources
- No repairs
- No loaners
- No service plans
- No accidental damage concerns
- No leak testing
- No caustic chemicals

**Maximize your cystoscopy volume with the aScope™ advantage:**

- **Single-use flexible cystoscope** – sterile and patient-ready out of the package
- Affordable, high performance, **all-in-one HD display** easily placed in or moved to any exam room
- Connects to standard fluid tubing and **works with a wide variety of common instruments**

aScope™ 4 Cysto | Ready When Needed




**Maximizing office-based cystoscopy with aScope 4 Cysto.**

Easy to implement for staff, high performance for physicians

- **Optimize schedule and workflow** to easily meet your evolving procedural and practice needs
- **Add-on cystoscopies** throughout the day without availability concerns
- Flexibility to **perform cystoscopies at new locations** previously limited by capital equipment or personnel
- **Reallocate staff resources** from time – consuming reprocessing to patient – focused care
- **Download images and videos** to share with patients or easily upload directly into your EMR to ensure data security
- **New cystoscope for every procedure** – no repairs, loaner scopes, or service plans ever required

**AMBU CONSULTANTS OFFER AT NO CHARGE:**

- Budget impact analysis reviewing current and future costs for cystoscopy
- Productivity improvement summary to evaluate how potential increases in procedures impact practice's revenue
- A single-use implementation assessment to review the approaches of utilization – from specific procedure rooms to full-office conversions





### aScope™ 4 Cysto

Single-Use, Digital, Flexible Cystoscope  
 Deflection 210/120 degrees  
 Insertion tube diameter 16.2 Fr  
 Instrument channel width 6.6 Fr  
 Depth of Field 3 - 100 mm  
 Field of View 120°

Product	Item No.	Quantity
aScope™ 4 Cysto, Standard Deflection	601001000	5 pieces
aScope™ 4 Cysto, Reverse Deflection	600001000	5 pieces

### aView™ 2 Advance

Portable, full-HD touchscreen monitor  
 Connectivity with PACS via DICOM (Wi-Fi / LAN)  
 Video and image capabilities  
 4+ hours of battery (or standard electr

Product	Item No.	Quantity
Ambu® aView™ 2 Advance	405011000	1 pieces



**A PATIENT-READY  
 CYSTOSCOPE AT A  
 MOMENT'S NOTICE.  
 NO CHEMICALS. NO  
 MACHINES. NO WAITING.**

The aScope 4 Cysto allows office-based urologists a powerful new tool for a wide array of cystoscopy procedures - both diagnostic and therapeutic – without concern for damage or downtime.

- Diagnostic (Urethra\Prostatic Urethra)
- Bladder Cancer Surveillance
- Ureteral Stent Removal
- Bladder Biopsies
- Bladder Fulguration
- Over-Active Bladder - Botox® Injection
- Urethral Stricture Dilation
- Difficult Catheter Placement



## THE LATEST INNOVATION from Ambu

A single-use solution from Ambu, the world leader in single-use flexible endoscopes, covering over 1 million procedures a year and growing.

A global organization with over 80 years of producing innovative medical devices and a mission to improve patient care, reduce costs, and enhance workflows.

**Produced nine flexible endoscopy systems over the last 10 years** and accelerating with new systems and specialties to help meet the demand for more safe and efficient health care - including our first ureteroscope in 2021/2022.

### Ambu's Vision of Partnership

Wide array of customer support programs including:

- On-site clinical training teams and continuing education courses
- Budget and health economic impact reports
- Risk management and compliance education
- Implementation support to maximize higher patient throughput

Solutions for...

- Pulmonology
- Critical Care
- Otolaryngology
- Urology
- Gastroenterology
- Anesthesiology

**Ambu: World Leader in  
Single-Use Endoscopy**

# Ambu

Ambu, Inc.  
6230 Old Dobbin Lane  
Columbia, MD 21045  
Tel. 800 262 8462  
Fax 800 262 8673  
ambuUSA.com

US: Rx only