

THE STERILE, SINGLE-USE DUODENOSCOPE SOLUTION

The Ambu Duodeno system consists of a sterile, single-use duodenoscope and the aBox Duodeno unit.

Remove the aScope Duodeno from its packaging, connect it to the aBox Duodeno and the system is ready.

The multipurpose system has an integrated rinsing function, and there is no need for an additional light source, saving you space and setup time.

The aBox Duodeno connects directly to a medical grade monitor and documentation systems in use.

It can be cleaned and disinfected with germicidal wipes.

Ambu® aBox™ Duodeno



KEY BENEFITS

Easy to setup



Seamless integration – the aBox Duodeno connects directly to a medical grade monitor and documentation systems in use

Video processor with integrated light control system and lens rinsing function

Multisource DVI and 2 X HD-SDI connections for video output



Ambu



SPECIFICATIONS

Power Supply	
Voltage	AC 120V / 230V
Frequency	50 / 60Hz
Consumption electric power	91 VA
Fuse Rating	2x 5A H 250V T

Size of the aBox Duodeno	
Dimensions	497 (D) x 380 (W) s 140 (H) mm
Weight	13 kg

ORDERING SPECIFICATIONS

aBox Duodeno Ordering Information	Quantity
485001000US	1
Used with:	
Ambu aScope Duodeno 482001000 Canada Only	Packs of 6
Ambu aScope Duodeno 482101000 US Only	Packs of 6



Ambu, Inc.
6721 Columbia Gateway Drive, Suite 200
Columbia, Maryland 21046
Tel. 800 262 8462
Fax 800 262 8673
ambuUSA.com

Questions? Contact your Ambu Representative
at 800-262-8462 or visit www.ambuUSA.com

