

Ambu® aScope™ Duodeno

STERILE. SINGLE-USE. PROBLEM SOLVED.

Sterile straight from the pack, Ambu aScope Duodeno is a single-use endoscope with side viewing optics, deflectable tip and an elevator to control the position of compatible endoscopic accessories.

With aScope Duodeno, you have a new scope with reliable performance every time, and there is no risk of patient cross-contamination.

What's more, there's no need for costly reprocessing or repair.



KEY BENEFITS

Sterile, single-use

No risk of patient cross-contamination

Cost effective



Eliminates the risk and cost of reprocessing

Familiar design and function to that of conventional duodenoscopes



Distal tip showing elevator in use with a sphincterotome

Ambu



SPECIFICATIONS - AMBU® ASCOPE™ DUODENO

Optical System	
Direction of view	Retro side viewing 6°
Field of view	130°
Insertion section	
Distal end outer diameter	φ13.7 mm
Insertion tube outer diameter	φ11.3 mm
Working length	1240 mm
Instrument working channel	
Inner diameter	φ4.2 mm (12.6 FR)
Bending section	
Angulation range	Up 120° / Down 90° Left 90° / Right 110°
Illumination method	Dual LEDs
Dimensions	
Weight	700 g

ORDERING SPECIFICATIONS

aScope Duodeno	Quantity
482001000	6
Used with Ambu aBox Duodeno - 485001000US	

Ambu

Ambu, Inc.
6230 Old Dobbin Lane
Columbia, MD 21045
Tel. 800 262 8462
Fax 800 262 8673
ambuUSA.com

CE 2797 US: Rx only