



Ambu

Ambu® Airway Management System

FROM INTUBATION TO ISOLATION

One system for all airways

ONE SOLUTION FOR ALL YOUR AIRWAY VISUALIZATION NEEDS

A connected visualization platform that aligns with Guidelines for the combination of video laryngoscopy and flexible bronchoscopy, while also offering continuous lung isolation monitoring for confident airway management.

Ambu® SureSight® Video Laryngoscope

Intubation - Maximize your focus on the patient



Stylets available in small, medium, and large

NEW! Ambu® SureSight® Mobile Video Laryngoscope



Ambu® aView™ 2 Advance

Endoscopy system - Your control center for key devices

Ambu® aScope™ 5 Broncho HD

Bronchoscopy - Next-level imaging and maneuverability



Ambu® aScope™ 4 Broncho

Bronchoscopy - Portability and essential performance



Ambu® VivaSight™ 2 DLT

OLV procedures - Built-in camera increases patient safety



Together with the visualization devices, Ambu laryngeal masks, face masks, resuscitators, and circuits comprise a comprehensive Ambu airway management solution.

The Ambu® video laryngoscopy solution

Ambu® aCart™ Compact

Everything within reach and ready to go – a compact, mobile workstation for the SureSight Connect video laryngoscope solution and your other key devices.

Ambu® aView™ 2 Advance
Ambu® aBox™ 2
(not shown here)

Plug-and-play system that serves as the control center for key Ambu airway visualization devices. Pair with an Ambu aCart Compact or mount directly to an IV pole to save space in a busy OR.

Ambu® SureSight® Video Laryngoscope

An intuitive, ergonomic and reliable video laryngoscope

- Audible click ensures that the blade is properly attached
- Ergonomic design for secure grip and optimal control
- Anti-fog lens for a clear and consistent image
- Easy-to-reach video/camera button



Comprehensive blade portfolio

VL blades for both routine and difficult intubations as well as emergent procedures

- Miller 0
- Miller 1
- MAC 3
- MAC 4
- Hyperangulated 1
- Hyperangulated 2
- Hyperangulated 2 Channeled
- Hyperangulated 3
- Hyperangulated 3 Channeled
- Hyperangulated 4

NEW! Ambu®
SureSight® Mobile
Video Laryngoscope

ONE SYSTEM FOR EVERY INTUBATION

The Ambu® SureSight® Video Laryngoscope is part of our trusted airway visualization platform. Having one uniform, intuitive, and familiar system helps you stay focused on the patient as you switch between devices and move through the airway management algorithm.

One system, multiple visualization devices

With the Dual View feature in the Ambu endoscopy system, you are able to simultaneously use two of our airway visualization devices without needing separate screens for difficult airways.

Streamlines work with patient imaging data

The Ambu endoscopy system offers connectivity with Electronic Health Records (EHR) to simplify the important task of sharing and transferring patient data.

[Learn more about the Ambu intubation solution](#)



WHY IT MATTERS:

Maintaining airway visualization throughout the procedure can help clinicians verify placement, identify displacement, and coordinate airway management.

The Ambu® bronchoscopy solution

The Ambu® aScope™ Broncho portfolio provides access to seven different bronchoscopes available in a variety of clinically appropriate sizes.

Ambu® aScope™ 5 Broncho

- 2.7/1.2 The world's first ultra-slim bronchoscope compatible with a 3.5 ETT that offers a working channel. Pair with SureSight Video Laryngoscopy for dual-view procedures, VivaSight DLT for optic assurance in OLV, and Fuji Uniblocker for simplicity in securing bronchial blocker placement.
- 4.2/2.2 A regular-sized working channel in a slim bronchoscope compatible with VivaSight™ DLT for continuous OLV visualization and dual-view with the SureSight® Video Laryngoscope. Also compatible with Fuji Uniblocker, simplifying bronchial blocker placements.
- 5.0/2.2 and 5.6/2.8 Larger outer diameter scopes with working channels designed for expedient airway clearance when it matters most.



Ambu® aScope™ 4 Broncho

- 3.7/1.2 (Slim)
- 5.0/2.2 (Regular)
- 5.8/2.6 (Large)



BRONCHOSCOPY BUILT FOR THE OR WORKFLOW

Trusted by surgical care teams around the globe for its performance, consistency and versatility.

aScope 4 Broncho: A complete portfolio for foundational bronchoscopy

Offers the essential performance and simplicity needed for foundational bronchoscopy procedures frequently performed in the OR and ICU.

aScope 5 Broncho: Access to high-performance advanced bronchoscopy starts here

Combines the capabilities for advanced diagnostic and therapeutic procedures with predictable scope performance for every case.

[Learn more about the Ambu bronchoscopy solutions](#)



WHY IT MATTERS:

Bronchoscopes are no longer just intubation tools in the OR. They support difficult airway management, DLT verification, airway inspections and secretion clearance – all within the same workflow

The Ambu® one-lung ventilation solution

Ambu® VivaSight™ 2 DLT

Built-in camera at the tip of the tracheal lumen enables continuous visualization during OLV procedures.



Fuji UNIBLOCKER™

A complete airway platform should support multiple OLV strategies – combining bronchial blockers with bronchoscopy gives clinicians additional flexibility when anatomy, procedure type, or airway complexity demands a tailored approach.



Ambu® aScope™ 4 and 5 Broncho

Comprehensive bronchoscope portfolio built to match the complexity and breadth of potential clinical scenarios during OLV procedures.



CONTINUOUS VISUALIZATION DURING ONE-LUNG VENTILATION

The Ambu® VivaSight™ 2 portfolio advances the safety and efficacy of One-Lung Ventilation (OLV) procedures with continuous, real-time visualization of tube placement. Unlike conventional double-lumen tubes (DLT) where tube positioning is blind without the use of a bronchoscope, VivaSight 2 features a built-in camera at the tip of the tracheal lumen, enabling clinicians to monitor positioning continuously throughout the procedure.

This continuous view allows immediate detection and correction of malpositioning—without compromising ventilation. In addition to clinical enhancements, workflow enhancements for bronchoscopy are provided when utilizing Dual View imaging by combining VivaSight 2 DLT with Ambu® aScope™ 4 or 5 Broncho, enabling simultaneous visualization of the DLT placement and bronchoscope on a single display.

For cases needing greater flexibility, such as complex airways or challenging anatomy, the Fuji UNIBLOCKER®, a bronchial blocker with a torque control shaft, can be used alongside a bronchoscope as an alternative lung isolation approach—expanding options within the Ambu airway visualization platform.

Because effective airway management isn't one-size-fits-all—Ambu provides you with the flexibility to choose the right approach for every patient.

Learn more about the Ambu VivaSight 2 OLV solution



WHY IT MATTERS:

Traditional OLV workflows often require clinicians to repeatedly stop and recheck tube positioning. Continuous visualization helps maintain awareness throughout the procedure while reducing workflow interruptions.

The Ambu® integrated display systems solution

Plug-and-play system that serves as the control center for key Ambu airway visualization devices



Ambu® aView™ 2 Advance Gen. 2

CAPTURE AND SHARE THE AIRWAY ENCOUNTER

Images and video from the Ambu platform connect directly for streamlined procedural records – especially critical in difficult airway cases.

[Learn more about the Ambu integrated solution](#)



WHY IT MATTERS:

Switching between disconnected airway systems can add complexity during critical moments. Integrated visualization platforms help simplify workflow by bringing airway imaging, documentation, and device coordination together.

ONE CONNECTED PLATFORM FOR THE ENTIRE AIRWAY WORKFLOW

In support of guidelines

From routine intubations to difficult airway management and one-lung ventilation, modern airway care depends on maintaining visualization, adaptability, and control at every stage of the procedure.

The Ambu airway visualization platform brings together video laryngoscopy, bronchoscopy, continuous lung isolation monitoring, and documentation within one connected workflow – helping clinicians move more seamlessly from intubation to inspection, reassessment, and ongoing airway management.

Because airway management isn't a series of disconnected tasks. It's one continuous clinical journey.

ORDERING SPECIFICATIONS

Ambu® SureSight® Video Laryngoscope

Item code	Product	Quantity
390511000	Ambu SureSight Connect US	1
390101000	Ambu SureSight Blade Miller Size 0 US	5
390111000	Ambu SureSight Blade Miller Size 1 US	5
390211000	Ambu SureSight Blade Hyperangulated Size 1 US	5
390221000	Ambu SureSight Blade Hyperangulated Size 2 US	5
390521000	Ambu SureSight Blade Hyperangulated Size 2 Channeled US	5
390231000	Ambu SureSight Blade Hyperangulated Size 3 US	10
390531000	Ambu SureSight Blade Hyperangulated Size 3 Channeled US	10
390241000	Ambu SureSight Blade Hyperangulated Size 4 US	10
390331000	Ambu SureSight Blade MAC Size 3 US	10
390341000	Ambu SureSight Blade MAC Size 4 US	10
321900010	Ambu SureSight Stylet Small	10
321900020	Ambu SureSight Stylet Medium	10
321900030	Ambu SureSight Stylet Large	10
390471000	SureSight Mobile Kit - US (SureSight Mobile w. Battery, Carry Case, Extra Battery, Battery Charger and Data Cable)	1
390481000	SureSight Mobile Mini Kit - US (SureSight Mobile w. Battery, Carry Case and Extra Battery)	1
390611000	SureSight Battery - US	1
390621000	SureSight Battery Charger - US	1
390631000	SureSight Data Cable - US	1
390641000	SureSight Carry Case - US	1

Ambu® aScope™ 5 Broncho / Ambu® aScope 5 Broncho™ Sampler Set

Item code	Product	Quantity
621001000US	Ambu® aScope™ 5 Broncho HD 5.0/2.2	5
622001000US	Ambu® aScope™ 5 Broncho HD 5.6/2.8	5
620001000US	Ambu® aScope™ 5 Broncho 4.2/2.2	5
624001000US	Ambu® aScope™ 5 Broncho 2.7/1.2	5
621002000US	Ambu® aScope™ 5 Broncho HD 5.0/2.2 Sampler Set	5
622002000US	Ambu® aScope™ 5 Broncho HD 5.6/2.8 Sampler Set	5
620002000US	Ambu® aScope™ 5 Broncho 4.2/2.2 Sampler Set	5

Ambu® aScope™ 4 Broncho / Ambu® aScope 4 Broncho™ Sampler Set

Item code	Product	Quantity
476101000	Ambu® aScope™ 4 Broncho Slim	5
477101000	Ambu® aScope™ 4 Broncho Regular	5
478101000	Ambu® aScope™ 4 Broncho Large	5
477501000	Ambu® aScope™ 4 Broncho Regular Sampler Set	5
478501000	Ambu® aScope™ 4 Broncho Large Sampler Set	5

Ambu® VivaSight™ 2

Item code	Product	Quantity
412351000	Ambu® VivaSight™ 2 DLT Kit 35 Fr - LEFT	1
412371000	Ambu® VivaSight™ 2 DLT Kit 37 Fr - LEFT	1
412391000	Ambu® VivaSight™ 2 DLT Kit 39 Fr - LEFT	1
412411000	Ambu® VivaSight™ 2 DLT Kit 41 Fr - LEFT	1

Ambu® Display Units and Carts

Item code	Product	Quantity
405011000US2	Ambu® aView™ 2 Advance Gen. 2 Displaying and Processing Unit	1
505001000	Ambu® aBox™ 2 Displaying and Processing Unit	1
KU.9713.800	Ambu® aCart™ Compact	1
KN.1215.991	Ambu® aCart™ Compact Vario Mount	1
KN.1245.800	Ambu® aCart™ Compact Universal Endoscope Holder	1



THE AMBU COMMITMENT

Ambu is on a mission to deliver airway management solutions that help you care for patients and save lives. As a single-use medical device company, we are also acutely aware of our obligation to pursue more sustainable solutions.

Taking leaps forward in sustainability

Advancing sustainability is one of the four strategic pillars of our company. We focus on circular design principles that reduce the carbon footprint of our products and on value-chain activities that help us move toward our net-zero emissions goal.

One concrete example is our pioneering use of bioplastic material – a mix of fossil-based and second-generation bio feedstock – in our endoscope handles, laryngeal mask cuff protector packaging, and other products.

By the end of 2024, we produced all our endoscopes with bioplastic material in the handles. In the future, we plan to introduce this material into other parts of our endoscopes. In addition, SureSight is phthalate-free and PVC-free.



Contact us for more information

Learn more about the Ambu airway visualization platform, including solutions for intubation, bronchoscopy and one-lung ventilation.

[Learn more about the Ambu Airway Management Platform](#)

Ambu, Inc.
6721 Columbia Gateway
DriveSuite 200
Columbia, MD 21046
Tel. 800 262 8462
Fax 800 262 8673
ambuUSA.com